**ENG
Maxxi Flow Sensor
User Instructions**

Intentions of Use:


The Maxxi Flow sensor is a respiratory sensor. It is placed under the patient's nostrils and produces a signal that is directly proportional to the temperature changes of air inhaled and exhaled during respiration. Maxxi Flow sensors generate a small analog electrical signal that provides a clear, reliable indication of respiration airflow. This sensor is intended to be used with polysomnography devices and adult patients. Intended to be use in a sleep laboratory, clinics, or hospitals.

Maxxi Flow Sensor Instructions for Use:

The respiration flow sensor is placed on the patient's upper lip, under the nostrils. The sensors are delicate and must be handled with care.

1. Gently bend the two nasal prongs forward to prevent them from entering the nostrils. Make a round bend slightly away from the sensor's body.
2. Align the tip of the oral prong with the center of the mouth/upper lip.
3. Draw the wires on either side of the sensor over the patient's ears and use the clasp under the chin to secure the wires in place.
4. Tape the wire on both cheeks to reduce strain and enhance reliability.
5. Check that clear, strong signals are being transmitted before starting the study.

Recommended polysomnography filter settings:


1. Low-cut filter: 0.1 Hz.
2. High-cut filter: 15 Hz.

Cleaning:

1. To clean the sensor and cables use with a cloth moisture with hospital grade water followed by cleaning/disinfect by wiping with a soft cloth dampened with hospital grade cleaner that is not corrosive to plastic or metal. The following materials may be used:
 - a. 70 % Isopropanol

- b. Super Sani-Cloth plus wipes
2. Do not over stretch the cable.
 3. Do not keep the connectors wet.
 4. Do not immerse in water.

Warnings and precautions:

1. Do not connect directly to the electric grid (RISK OF ELECTRIC SHOCK).
 2. THIS DEVICE IS NOT TO BE USED IN MR ENVIRONMENT.
 3. Caution: Federal law restricts this device to sale by or on the order of a physician.
 4. Medical product NON-STERILE. - Do not use the product if the package is previously opened or violated.
 5. Maxxi Flow sensor is for medical professional use only such, physicians and certified sleep technologists (BRPT or CPSGT).
 6. Maxxi Flow sensor may be used only in conjunction with FDA cleared recording device or approved by a similar regulatory authority in your country.
- 
7.  For correct recycling process, please follow local regulations to dispose of the electrical material, if applicable.
 8. Discontinue use if sensor shows signs of wear, damage, or exposed metal (example: plastic cracks, missing connector, missing thermocouple element, debris or corrosion, loose connector, non-factory color).
 9. During eventual electromagnetic disturbances, which affect the essential performance of the device, there will be a noticeable degradation of the physiological signals displayed during registration, this degradation differs from the physiological signals and should be interpreted / considered as artifacts (interferences), in view of this situation we recommend that the user disregard the record until the disturbances is ceased.

Operating/Storage Conditions:

Use type	Multiple patient multiple use
Operating conditions	15°C (50°F) - 30°C (95°F)
Storage temperature	-20°C (-4°F) - 60°C (140°F)
Operating and storage humidity	5% - 95% (non-condensing)

If any physical disturbance and/or chemical breakdown happens during the storage period, do not use the product.

Warranty Period:

Neurovirtual USA warrants the Interface module to be free of defects in materials and workmanship for a period of 6 months from the date purchased.

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Neurovirtual USA Inc.