

Now you can count on the best online store to provide all of your PSG/EEG supplies.

www.sleepvirtual.com/products

Use the "welcome" coupon code in your first purchase to get a 10% discount off your entire order.

Thermocouple Flow Sensor

The Sleepvirtual Thermocouple is a durable sensor with stiff bendable prongs, which allow you to adjust the sensor for optimal patient comfort. Comprised of three heat-sensitive elements, the Thermocouple Flow Sensor produces a signal that is directly proportional to the temperature changes of air inhaled and exhaled during respiration.



RIP/Inductive Technology

Sleepvirtual RIP/Inductive Technology has one of the lowest setup costs available. Sleepvirtual bands are available in neonatal to XXL blue or white. They are comfortable for the patient and provide high-quality signals with little-to-no adjustment throughout the night.



Pressure Transducer

The Sleepvirtual Pressure + Snore Sensor employs a piezo membrane to generate a small voltage in response to pressure changes and acoustical input.



Snore Sensors

Using a Piezo-Electric Crystal, the SleepSense Piezo Snore Sensor converts vibrations in the neck caused by snoring to a small analog voltage that provides a clear, reliable indication of the presence of snoring.



Body Position Sensors

Sleepvirtual Body Position Sensors provide the most accurate report of body position. Small and robust, they will deliver years of high-quality service. They work with most available PSG systems.



Pediatric Sensors

Sleepvirtual has a wide variety of sensors specifically designed for infants and pediatrics. They are comfortable for pediatric patients, easy to apply and produce high-quality output signals.



Flex Sensors SpO2

Nonin Medical's Flex Sensor System combines the convenience and performance of a disposable sensor with the economics of a reusable sensor.



Now you can count on the best online store to provide all of your PSG/EEG supplies.

www.neurovirtual.com/products

Use the "welcome" coupon code in your first purchase to get 10% of discount in your entire order.

Conductive Paste

Ten20 contains the right balance of adhesiveness and conductivity, enabling the electrodes to remain in place while allowing the transmittance of electrical signals.



Skin Prep Gel

Nuprep Skin Prep Gel effectively lowers impedance to improve your tracings. Its mildly-abrasive formula improves conductivity and helps achieve maximum efficiency with your equipment.



EEG Cup Electrodes

Used for electroencephalographic testing, the gold-plated electrode—reusable with center hole—allows for easy fill of electrode paste after adhesion of cup to the scalp. Available in 48", 60", 96" and 120" lead wires length.

Other Application Items



EEG Templates



Electrocap



Cannulas



Snap Electrodes



Button Snap Electrodes



Head Net



Adhesive Tapes



Markers