



## Product Specifications BWFlash



- Meets the most rigorous methodologies for photo stimulation (0.7 joules of intensity at 30 centimeters from the patient's face).
- Intense luminosity, thanks to multiple high intensity and efficient LEDs, without heat.
- Independent control on the panel, without the need to connect to a Notebook.
- 10 cm diameter diffuser for uniform stimulation.
- Programmable through BWFlash software (compatible with Windows® 10).
- Integration available both to the Neurovirtual electroencephalograph line or other devices.

## Technical Specifications BWFlash



Characteristics	Models / Values
Frequency of stimulation	1 - 60 Hz
Power of stimulation	10 levels - up to 13000 lux peak at 6 inches
Models of operation	Stand alone mode
	Software controlled mode
	Trigger out mode
Power supply	5VDC
Connection Type	USB Type B
Trigger out connection	P2 3.5 mm

Dimensions BWFlash			
Length	Height	Width	Weight
152mm (6")	127mm (5")	114mm (4.5")	370 gr (13oz)